

# ONYX ALL-IN-ONE (AIO) MEDICAL PC

## Cleaning Instructions

### How to clean the entire AIO PC for infectious control

- Step 1. Use a disposable alcohol based anti-bacterial wipe (approx. 70% alcohol)
- Step 2. Begin by cleaning screen and OSD keys at the front of the AIO PC
- Step 3. Clean the remainder of the rear of the AIO PC.

### How to clean the screen of AIO PC for general house keeping

- Step 1. Use either a microfibre cloth with an alcohol based cleaning spray (approx. 70% alcohol) or disposable alcohol based cloth to do general housekeeping on the screen.

### Onyx Chemical Cleaning Guide

Please download the Onyx Chemical Cleaning Guide from [micromaxhealth.com](http://micromaxhealth.com) for a full list of acceptable chemicals that can be used on Onyx devices.