

IPTEC Level 2 Surgical Masks

IPTEC masks are manufactured by Sys-mac Automation Engineering in Singapore with dedication to quality. Manufactured under strict compliance of ISO 13485:2016 they provide a **high quality** and **value for money** product.

This product is registered with the TGA. Micromax Pty Ltd is the proud sponsor. **ARTG: 341086.**

PRODUCT EXAMPLES:



50 masks per box, vacuum sealed to prevent microbial growth. Single-Use product only. Non-sterile.



SURGICAL MASK FEATURES:

- Bacterial Filtration Efficiency (BFE) 99%
- High Breathability with excellent differential pressure value
- Fluid resistant
- Flammability tested
- Non-woven & Hypoallergenic
- Soft Ear Loop
- Nose-bar Adjustable for perfect fit
- 3 Layers of protection.



TECHNICAL DETAILS:

- 3-Ply pleated design, Ear Loop
- Single-Use, Non-sterile Face mask
- AS 4381: 2015 Level 2
- ARTG ID: 341086
- EN14683 Type IIR.



To discuss your specific requirements with our team or for a free, no obligation quote contact us on **1300 760 699** or email **info@micromaxhealth.com.au**

