

# The Benefits on Using XXL Battery

*Safest, Most elegant & Cost Effective  
Mobile Powered Solution*

Items	XXL Battery	XL Battery
Powerful	<b>+40%</b>	
Runtime (50% video,50% idle)	<b>18 Hours</b>	12 Hours
Longevity	<b>5 Years</b>	3 Years

- Improved Nursing Work-flow

- *No handling of batteries*
- *No sanitising of batteries*
- *Prevent battery dropping & associated hazards with broken casings*

- Similar Performance to Expensive Powered Carts

- *Trickle feed power for onboard charging just like the mobile devices*
- *Simply plug into mains power when required for charging batteries*



- **2 XXL Battery Kit. Li-ion. 10.8V. 12060mAh**
- **2 Years 1000 cycle warranty**
- **130Wh battery capacity**
- **Smart LED: battery status indicator for end user**