




## Smart View Medical Station

# ZEUS-228

22" Fanless Smart View 6th Generation Core i7  
Multi-Touch Medical Station



## Features


- Intel® 6th Generation Core™ i7-6600U 3.4GHz Dual Core Processor
- 22" 250nits full HD high contrast MVA LCD
- Support DDR4 SODIMM up to 32GB
- Capacitive multi-touch screen/5 wire resistive touch screen
- Easy cleaning flat screen
- 10 programmable smart function keys with dual colours LED status indication
- LAN /USB (optional)/COM(optional) medical isolation
- Super speed 10Gb/s USB 3.1 Port
- PCI Express[x4] x 1, PCI Express[x1] , PCI x1(optional)
- Trusted platform module/smart card reader(optional)/RFID reader(optional)
- Remote management: Intel active management technology
- +12V to +24V wide range DC input
- 22Wh backup battery (optional)

## Specifications

Main Specifications	
Processor	Intel® 6th Generation Core i7-6600U 3.4GHz Dual Core
System Memory	DDR4 SODIMM up to 32GB
Expansion	PCI Express[x4] x 1, PCI Express[x1] , PCI x1(optional)
Storage Disk Drive	2.5" SATA Hard Disk Drive/SSD x 1 Slim DVD-RW Drive x 1(optional) M.2 SSD x 1
Security	Trusted Platform Module, Smart Card reader(optional), RFID(optional)
Wireless Communication	802.11a/b/g/n/ac(optional), Bluetooth 4.0 (optional)
Speaker	5W x 2
Function Key	Speaker Volume Up/Down, LCD Brightness Up/Down, LCD On/Off, Touch Screen On/Off, DICOM Mode(optional), Reading Light On/Off
Power Requirement	DC12~24V, 22Wh back up battery(optional)
OS Support	Windows® 7, Windows® 8.1, Windows® 10, Linux®
Display	
Size	22" LCD
Resolution	1920 x 1080
Luminance	250 nits
View Angle	178°(H)/178°(V)
Contrast Ratio	3000:1
Back Light Life Time	30,000 hours
Touch Screen	Capacitive multi-touch/5-wire resistive
I/O	
USB	USB 3.0 x 4, USB 3.1 x 2(side I/O) Isolated USB 2.0 (optional)
Serial Port	Isolated RS-232 x 2(optional)
Ethernet	Isolated gigabit LAN x 2
Video Out	Display port x 2 (optional)
Mechanical and Environmental	
Power Consumption	Full loading: 55 Watts
Operating Temperature	0°C ~ 40°C(32°F ~ 104°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Mounting	VESA 75/100 mm
Degree of Protection	IP65 in the front, IPX1
Net Weight	8kg (17.6lb)
Package Dimension (W x D x H)	750 x 220 x 490mm
Certifications	CE, FCC Class B Passed 18,cULus Mark/UL+VTC Report+ISO 14971(CAN/CSA-C22.2NOo.60601-1-08; ANSI/AAMI ES60601-1:2005, EC 60601-1-1:2005, EN 60601-1:2006 ISO 14971)

## 22" Fanless Smart View 6th Generation Core i7 Multi-Touch Medical Station

### Dimension / Unit: mm



### Ordering Information

- |  |   |
|--|---|
| » ZEUS-228ST-A1-1010<br>Medical Station.i7-6600U.4GB RAM.22" LCD.12~24V.Capacitive Touch | » ZEUS-228ET-A1-1010<br>Medical Station.i7-6600U.4GB RAM.22" LCD.12~24V.Resistive Touch |
| » ZEUS-228ST-A2-1010<br>Medical Station.i5-6300U.4GB RAM.22" LCD.12~24V.Capacitive Touch | » ZEUS-228ET-A2-1010<br>Medical Station.i5-6300U.4GB RAM.22" LCD.12~24V.Resistive Touch |
| » ZEUS-228ST-A3-1010<br>Medical Station.i3-6100U.4GB RAM.22" LCD.12~24V.Capacitive Touch | » ZEUS-228ET-A3-1010<br>Medical Station.i3-6100U.4GB RAM.22" LCD.12~24V.Resistive Touch |

### Optional Accessories

- |   |  |
|---|--|
| » 1255301202<br>Medical Power Adapter, 120W, 24V                        | » OPM-C02C-A3<br>One Isolated USB 2.0 + Two Isolated R-232 PCIe Card |
| » OPM-C06W-A6<br>Mini Card WLAN Kit, IEEE 802.11a/b/g/n + Bluetooth 4.0 | » OPM-V03B-A1<br>One Display Port + Two RS-232 Kit                   |
| » OPM-T003-A4<br>Slim Slot-in DVD-RW Drive                              | » OPM-V03B-A2<br>One Display Port Kit                                |
| » OPM-T006-A1<br>Smart Card Reader Kit                                  | » OPM-T012-A2<br>Two RS-232 Kit                                      |
| » OPM-P03C-A4<br>22Wh Backup Battery                                    | » OPM-V04B-A2<br>DICOM Kit   |
| » OPM-C02C-A4<br>Two Isolated RS-232 PCIe Card                          | » OPM-R02X-A1<br>PCI Express [X1] Riser Kit                          |
|   | » OPM-R01P-A1<br>PCI Riser Kit                                       |